

## Model 30<sup>™</sup> Pneumatonometer

Precise tonometry, ocular pulse amplitude, and tonography measurements.





## Model 30<sup>™</sup> Pneumatonometer

The Reichert<sup>®</sup> **Model 30<sup>™</sup> Pneumatonometer** is an easy to use instrument that provides fast and accurate tonometry and optional tonography functions. The probe tip, which floats on an air bearing, is gently touched to the anesthetized cornea and a precisely regulated flow of filtered air applies force to the tip. A small (5 mm diameter) fenestrated membrane permits the air to flow through vents in the tip until it conforms to the shape of the cornea. Increasing air pressure is continually applied to the cornea until the force being applied is equal to the pressure in the anterior chamber. When these forces are in balance a pneumatic sensor records the intraocular pressure. In pulse tonometry mode, the ocular pulse waveform is charted and recorded along with IOP.

### **Unique Methodology**

The Model 30 Pneumatonometer measures intraocular pressure noninvasively through applanation tonometry. A gentle, floating pneumatic sensor touches the surface of the anesthetized cornea with the exact amount of applanation force required to take the measurement.

### Easy to Use

Doctors and technicians agree that the Model 30 Pneumatonometer is easy to use. The touch screen interface features an easy-to-use icon-based operating system enabling even inexperienced users to make measurements in seconds. The full color display guides the operator through the measurement process and displays readings. The optional footswitch (available with the tonography kit – catalog number 16033) makes obtaining measurements even easier, freeing the operator to concentrate on handling the sensor probe. The optional Model 30 Pneumatonometer PC Software enables data to easily be transferred to a PC.

### **Multi-Modality**

The Model 30 Pneumatonometer provides ocular pulse amplitude measurements in addition to tonometry measurements. The optional tonography mode enables measurement of the aqueous outflow efficiency of the trabecular meshwork.

### **Clinically Trusted**

The Model 30 Pneumatonometer is trusted worldwide and is frequently the instrument of choice of specialists and researchers. Research indicates that measurements obtained with the Model 30 Pneumatonometer are less affected by corneal properties than other methods of tonometry.



## **Specifications**

Dimensions	14"W x 10.5"D x 5.25"H (35.6 x 26.7 x 13.3 cm
Weight	8.5 lbs (3.83 kg)
Air Pump	14 PSI, 3,000 hour life
Probe	Pneumatic
IOP Range	5 - 80 mmHg
Display	4.3 inch, touch screen, color LCD TFT
Power	100 V, 120 V or 230 V

#### **Catalog Numbers**

16030	Model 30 Pneumatonometer 120 V
16031	Model 30 Pneumatonometer 230 V
16032	Model 30 Pneumatonometer 100 V
16033	Tonography Kit: includes 10 g weight and footswitch

#### Accessories\*

	ichert com to purchase opline
232373	Calibration Verification Kit
232346	Filter Maintenance Kit
230678	Tip and Membrane Assembly 10/pk
230677	Tip and Membrane Assembly 3/pk
230676	Tip and Membrane Assembly 1/pk
232349	Probe Assembly (Includes Tip and Membrane assembly)

\*Visit reichert.com to purchase online



Proudly designed & assembled in the USA.

# For more information visit us online at www.reichert.com or call 1-888-849-8955.

#### Important Notice

Federal law restricts this device to sale by or on the order of a physician. Model30-1021 @2021 Reichert, Inc. Specifications subject to change without notice. Reichert and Reichert Technologies are registered trademarks of Reichert, Inc. Model 30 is a trademark of Reichert, Inc. AMETEK is a registered trademark of AMETEK, Inc. Reichert, Inc. 3362 Walden Avenue, Suite 100 Depew, NY 14043, USA Toll Free: 1-888-849-8955 • Tel: +1 716-686-4500 www.reichert.com • Iso 13485 Certified Reichert products are designed and manufactured under quality processes meeting ISO 13485 requirements.