



IMPROVED EFFICIENCY

ONE STEP MVR CANNULA SYSTEM 23, 25 & 27 GAUGE







ANOTHER FIRST IN IMPROVING VR SURGERY FROM DORC

- RÉVOLUTIONARY NEW MVR BLADE. REQUIRES MORE THAN 50% LESS INSERTION PRESSURE.
- UNIQUE TRANSLUCENT, NEWLY DEVELOPED CLOSURE VALVE PROVIDES ENHANCED VISIBILITY WHEN ENTERING INSTRUMENTS.

MORE THAN 50% LESS INSERTION PRESSURE REQUIRED







NEW CLOSURE VALVE PROVIDES ENHANCED VISIBILITY WHEN ENTERING INSTRUMENTS












CONTINUUM RANGE, THE 23/25/27G TDC CUTTER

A NEW GENERATION **OF VITRECTOMY CUTTERS FOR QUICK** AND SAFE REMOVAL **OF VITREOUS**











UNIQUE FEATURES AND BENEFITS:


- SAFER: DOUBLE CUTTING SPEED REDUCES TRACTION ON SURROUNDING TISSUE.


 MORE EFFICIENT: DUTY CYCLE OF 92% CREATES CONSTANT ASPIRATION FLOW INDEPENDENT FROM CUTTING SPEED.

 INCREASED ASPIRATION FLOW UP TO THREE TIMES MORE THAN NORMAL RATE.

 BETTER: REDUCED 'SURGE' TURBULENCE AT ASPIRATION PORT FROM CONSTANT OPENING OF THE ASPIRATION PORT.

 STRONGER: CUTTING PERFORMANCE FOR REMOVAL OF STRONG PVR MEMBRANES / LENS FRAGMENTS, DUE TO 42 PSI WORKING PRESSURE.








HIGH MYOPIC INSTRUMENTS


ENLARGING THE OPTIONS FOR TREATMENT OF HIGH MYOPIA




HIGH MYOPIC INSTRUMENTS

- FORCEPS AND BACKFLUSH INSTRUMENTS ARE ESSENTIAL FOR OPERATING ON THE HIGH MYOPIC EYE, AND YET, STANDARD INSTRUMENTS ARE OFTEN TOO SHORT TO REACH THE POSTERIOR POLE.
- THAT'S WHY DORC HAS DEVELOPED A UNIQUE RANGE OF INSTRUMENTS FOR TREATING HIGH MYOPIC EYES.

HIGH MYOPIC FORCEPS




HIGH MYOPIC BACKFLUSH INSTRUMENTS





23/25 GAUGE INSTRUMENTS & ACCESSORIES

25G CURVED STEPPED LASERPROBES



25G RETRACTABLE SWEEPER







CLAUS ECKARDT, FRANKFURT (GERMANY)

23/25 GAUGE INSTRUMENTS & ACCESSORIES





ILLUMINATION PROBE WITH

PHOTON CONNECTOR, INCLU-

DING ILLUMINATED SCLERAL

DESPRESSOR (23G/0.6 mm)



SHIFLDED TOTALVIEW ENDO-ILLUMINATION PROBE WITH PHOTON CONNECTOR, INCLU-DING ILLUMINATED SCLERAL DESPRESSOR (25G/0.5 mm)





NEW DISPOSABLE I/A HANDPIECES

DISPOSABLE COAXIAL HANDPIECES









I/A HANDPIECE, 45° ANGLED TIP FOR 2.0 mm INCISION, POLISHED ASPIRATION HEAD (BOX/10, STERILE)

DISPOSABLE BIMANUAL I/A SET



- BIMANUAL I/A SET: (BOX/10, STERILE)

 ASPIRATION HANDPIECE 0.8 mm, 21G WITH TEXTURED TIP &

 IRRIGATION HANDPIECE 0.8 mm, 21G WITH DOUBLE SIDED IRRIGATION PORTS

 FOR 1.2 mm INCISION



Scheijdelveweg 2, 3214 VN Zuidland The Netherlands

Phone: +31 181 45 80 80 Fax: +31 181 45 80 90 E-mail: sales@dorc.eu

Dutch Ophthalmic USA

Phone: +1 800 75 DUTCH or +1 603 778 69 29

sales@dutchophthalmicusa.com

D.O.R.C. Latin America Phone: +1 305 663 56 62 E-mail: la-sales@dorc.eu

D.O.R.C. Middle East, Africa and India

Phone: +971 509 11 65 87 E-mail: cvds@dorc.eu

D.O.R.C. Scandinavia AB

Phone: +46 19 12 15 14 E-mail: order@dorc.eu

D.O.R.C. France S.A.R.L. Phone: +33 1 53 58 14 73 E-mail: ventes@dorc.eu

D.O.R.C. Deutschland GmbH +49 211 52 29 15 80 E-mail: bestellung@dorc.eu

D.O.R.C. Austria GmbH +43 1 59 99 96 30 E-mail: info-at@dorc.eu

D.O.R.C. UK Limited

Phone: +44 80 00 51 83 77 E-mail: orders-uk@dorc.eu

